



ClearFit[®]


PATIENT SPECIFIC CRANIAL IMPLANTS



In a bedside and outpatient setting ClearFit implants allow for post operative imaging using ultrasound.

l o n g @ v i t i


Standard ClearFit



Single Stage ClearFit




Complex ClearFit



PROCESS

- 1 **Patient Scan**- Upload the patient's CT scans and order form
- 2 **Interactive Planning**- You'll be invited to an interactive planning session with our engineers to discuss your preferences and plans for the case
- 3 **Design Approval**- Design proposals are sent via PDF for surgeon review and approval
- 4 **Production**- Manufacturing begins at 3D printing and continues through sterilization
- 5 **Delivery**- Sterile implants are delivered prior to scheduled surgery date



www.longevity.com
sales@longevity.com



Patented Technology